


# EZO Complete-pH™

USB pH meter

**Datasheet for engineers**

**ISO 10523 Compliant**

(determination of pH)



Reads	<b>pH</b>	Calibration	<b>1, 2, 3 point</b>
Normal range	<b>.001 – 14.000</b>	Recalibration frequency	<b>~8 – 12 months</b>
Extended range	<b>-1.600 – 15.600</b>	Data protocol	<b>Serial data through FTDI virtual comport</b>
Accuracy	<b>+/- 0.002</b>	Temp compensation	<b>Yes</b>
pH reading time	<b>800ms</b>	Data format	<b>ASCII</b>
Supported probes	<b>Any type &amp; brand</b>	Ingress protection	<b>IP62</b>




Written by Jordan Press  
Designed by Noah Press

**PATENT PROTECTED**

This is an evolving document, check back for updates.

The EZO Complete-pH™ has all the features of this bench top meter.



- 1 Three decimal pH reading
- 2 Temperature used for reading
- 3 Calibration slope
- 4 Extended range capability
- 5 Immediate reading
- 6 Timed readings
- 7 Set device name
- 8 Voltage usage
- 9 Multi-point variable calibration
- 10 Manual Temperature compensation


The EZO Complete-pH™ is compatible with any brand of pH probe.

# Table of contents


Power consumption	4
Absolute max ratings	4
Interference free	5
Ingress protection – IP62	5
Default state	7

LED color definition	8
Receiving data from device	9
Sending commands to device	10
Simple serial monitor	11
<b>UART quick command page</b>	<b>12</b>
LED control	13
Find	14
Continuous reading mode	15
Single reading mode	16
Calibration	17
Export calibration	18
Import calibration	19
Slope	20
Extended pH scale	21
Temperature compensation	22
Naming device	23
Device information	24
Response codes	25
Reading device status	26
Sleep mode/low power	27
Factory reset	28

Calibration theory	29
Understanding pH Slope	34
Datasheet change log	37
Firmware updates	38
Warranty	39



# The EZO Complete-pH™ consists of 3 major components.



## Power consumption


5V USB	MAX	STANDBY	SLEEP
	37.0 mA	36.8 mA	22.6 mA

## Absolute max ratings

Parameter	MIN	TYP	MAX
Storage temperature	-65 °C		125 °C
Operational temperature	-40 °C	25 °C	85 °C

# Interference free


The EZO complete readings are unaffected by other sensors in the same water.



# Ingress protection – IP62


The EZO Complete-pH™ is dust proof and resistant to splashing water.

Two areas of concern are the **USB connector** and the **probe connector**.




# Ingress protection – IP62


An electrical short can occur if water enters the USB connector. A USB short could permanently damage the EZO-Complete. A USB short is not covered under warranty.



A connector short can occur if water enters the SMA connector. A connector short will cause the pH readings to pin to 0, 14, or the probe will respond slowly to changes in pH. A connector short is reversible and will not damage the EZO-Complete. However, frequent shorts will eventually damage the pH probe.




The SMA connector is part of your probe; Nothing should be in contact with this part.




**Do not use this  
as a mounting point!**

# Default state


**Baud** 9,600  
**Readings** continuous  
**Speed** 1 reading per second




# LED color definition




**Green**  
UART standby




**Cyan**  
Taking reading



**Purple**  
Changing  
baud rate



**Red**  
Command  
not understood



**White**  
Find

**5V**

LED ON  
**+2.2 mA**

**3.3V**

**+0.6 mA**

## Settings that are retained if power is cut

- Calibration
- Continuous mode
- Device name
- Enable/disable response codes
- LED control
- Protocol lock

## Settings that are **NOT** retained if power is cut

- Find
- Sleep mode
- Temperature compensation



# Receiving data from device


2 parts

ASCII data string

Command

Carriage return <cr>

Terminator



## Advanced

ASCII: 9 . 5 6 0 <cr>

Hex: 39 2E 35 36 30 0D

Dec: 57 46 53 54 48 13

# Sending commands to device


2 parts

**Command (not case sensitive)**

ASCII data string

**Carriage return <cr>**

Terminator



(COM Port numbers are determined by the computer)

## Advanced

ASCII: **S** **I** **e** **e** **p** **<cr>**


Hex: **53** **6C** **65** **65** **70** **0D**

Dec: **83** **108** **101** **101** **112** **13**

# Looking for a simple serial monitor for debugging?

**Termite:** a simple RS232 terminal


[Click here to download](#)



## Settings

**Baud** 9600      **Parity** none  
**Data bits** 8      **flow control** none  
**Stop bits** 1      **Forward** none  
**Transmitted text** append CR

Enter commands here



Command responses are shown here

# Command quick reference

All commands are ASCII strings or single ASCII characters.

Command	Function		Default state
<b>C</b>	enable/disable continuous reading	pg. 15	enabled
<b>Cal</b>	performs calibration	pg. 17	n/a
<b>Export</b>	export calibration	pg. 18	n/a
<b>Factory</b>	enable factory reset	pg. 28	n/a
<b>Find</b>	finds device with blinking white LED	pg. 14	n/a
<b>i</b>	device information	pg. 24	n/a
<b>Import</b>	import calibration	pg. 19	n/a
<b>L</b>	enable/disable LED	pg. 13	enabled
<b>Name</b>	set/show name of device	pg. 23	not set
<b>pHext</b>	enable/disable extended pH scale	pg. 21	disabled
<b>R</b>	returns a single reading	pg. 16	n/a
<b>Sleep</b>	enter sleep mode/low power	pg. 27	n/a
<b>Slope</b>	returns the slope of the pH probe	pg. 20	n/a
<b>Status</b>	retrieve status information	pg. 26	enable
<b>T</b>	temperature compensation	pg. 22	25°C
<b>*OK</b>	enable/disable response codes	pg. 25	enable

# LED control

## Command syntax

L,1 <cr> LED on **default**

L,0 <cr> LED off

L,? <cr> LED state on/off?

## Example

## Response

L,1 <cr>


\*OK <cr>

L,0 <cr>


\*OK <cr>

L,? <cr>

?L,1 <cr> or ?L,0 <cr>  
\*OK <cr>



L,1



L,0

# Find

## Command syntax

This command will disable continuous mode  
Send any character or command to terminate find.


Find <cr> LED rapidly blinks white, used to help find device

## Example

## Response

Find <cr>

\*OK <cr>



# Continuous reading mode

## Command syntax

- C,1 <cr>** enable continuous readings once per second **default**
- C,n <cr>** continuous readings every n seconds (n = 2 to 99 sec)
- C,0 <cr>** disable continuous readings
- C,? <cr>** continuous reading mode on/off?

## Example

## Response

**C,1 <cr>**

**\*OK <cr>**  
**pH (1 sec) <cr>**  
**pH (2 sec) <cr>**  
**pH (n sec) <cr>**

**C,30 <cr>**

**\*OK <cr>**  
**pH (30 sec) <cr>**  
**pH (60 sec) <cr>**  
**pH (90 sec) <cr>**

**C,0 <cr>**

**\*OK <cr>**

**C,? <cr>**

**?C,1 <cr> or ?C,0 <cr> or ?C,30 <cr>**  
**\*OK <cr>**

# Single reading mode

## Command syntax

A single reading takes 800ms


**R** <cr> takes single reading

## Example


**R** <cr>

## Response


**9.560** <cr>  
**\*OK** <cr>



**Green**  
**Standby**



**Cyan**  
**Taking reading**



**Transmitting**



**800 ms**



# Calibration

## Command syntax

Issuing the cal,mid command after the EZO™ pH circuit has been calibrated, will clear the other calibration points. Full calibration will have to be redone.

**Cal,mid,n** <cr> single point calibration at midpoint  
**Cal,low,n** <cr> two point calibration at lowpoint  
**Cal,high,n** <cr> three point calibration at highpoint  
**Cal,clear** <cr> delete calibration data  
**Cal,?** <cr> device calibrated?

## Example

## Response

**Cal,mid,7.00** <cr>

**\*OK** <cr>

**Cal,low,4.00** <cr>

**\*OK** <cr>

**Cal,high,10.00** <cr>

**\*OK** <cr>

**Cal,clear** <cr>

**\*OK** <cr>

**Cal,?** <cr>

**?Cal,0** <cr> or **?Cal,1** <cr> or  
one point  
**?Cal,2** <cr> or **?Cal,3** <cr>  
two point three point  
**\*OK** <cr>

# Export calibration

## Command syntax

Export: Use this command to download calibration settings

Export,? <cr> calibration string info

Export <cr> export calibration string from calibrated device

## Example

## Response

Export,? <cr>

10,120 <cr>

### Response breakdown

10, 120

# of strings to export

# of bytes to export

Export strings can be up to 12 characters long, and is always followed by <cr>

Export <cr>

59 6F 75 20 61 72 <cr> (1 of 10)

Export <cr>

65 20 61 20 63 6F <cr> (2 of 10)

(7 more)

⋮

Export <cr>


6F 6C 20 67 75 79 <cr> (10 of 10)

Export <cr>

**\*DONE**

Disabling \*OK simplifies this process

Export <cr>



# Import calibration

## Command syntax

Import: Use this command to upload calibration settings to one or more devices.

Import,n <cr> import calibration string to new device

## Example

Import, 59 6F 75 20 61 72 <cr> (1 of 10)

Import, 65 20 61 20 63 6F <cr> (2 of 10)

⋮

Import, 6F 6C 20 67 75 79 <cr> (10 of 10)

## Response


\*OK <cr>

\*OK <cr>


⋮

\*OK <cr>

Import,n <cr>



\*OK <cr>  
system will reboot



\* If one of the imported strings is not correctly entered, the device will not accept the import, respond with \*ER and reboot.

# Slope

## Command syntax

After calibrating a pH probe issuing the slope command will show how closely (in percentage) the calibrated pH probe is working compared to the "ideal" pH probe.

`Slope,?` `<cr>` returns the slope of the pH probe

## Example

```
Slope,? <cr>
```

## Response

```
?Slope,99.7,100.3, -0.89 <cr>  
*OK <cr>
```

## Response breakdown

`?Slope,`

**99.7**


↑  
99.7% is how closely the slope of the **acid** calibration line matched the "ideal" pH probe.

**100.3**

↑  
100.3% is how closely the slope of the **base** calibration matches the "ideal" pH probe.

**-0.89**

↑  
This is how many millivolts the zero point is off from true 0.



# Extended pH scale

Very strong acids and bases can exceed the traditional pH scale. This command extends the pH scale to show below 0 and above 14.

Lowest possible reading: **-1.6**  
Highest possible reading: **15.6**

## Command syntax

- `pHext,0 <cr>` extended pH scale off (0–14) **default**
- `pHext,1 <cr>` extended pH scale on (-1.6–15.6)
- `pHext,? <cr>` extended pH scale on/off?

## Example

## Response

`pHext,1 <cr>`


`*OK <cr>`

`pHext,0 <cr>`

`*OK <cr>`

`pHext,? <cr>`

`?pHext,1 <cr>` or `?pHext,0 <cr>`



# Temperature compensation

## Command syntax

Default temperature = 25°C  
Temperature is always in Celsius  
Temperature is not retained if power is cut

**T,n** <cr> n = any value; floating point or int

**T,?** <cr> compensated temperature value?

**RT,n** <cr> set temperature compensation and take a reading

## Example

## Response

**T,19.5** <cr>


**\*OK** <cr>

**RT,19.5** <cr>


**\*OK** <cr>  
**8.91** <cr>

**T,?** <cr>

**?T,19.5** <cr>  
**\*OK** <cr>



→  
**T,19.5** <cr>



# Naming device

## Command syntax

Do not use spaces in the name

Name,n <cr> set name

Name, <cr> clears name

Name,? <cr> show name

n =

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16

Up to 16 ASCII characters

## Example

## Response

Name, <cr>

\*OK <cr> name has been cleared


Name,zzt <cr>

\*OK <cr>

Name,? <cr>


?Name,zzt <cr>  
\*OK <cr>

Name,zzt



\*OK <cr>

Name,?



?Name,zzt <cr>  
\*OK <cr>

# Device information

## Command syntax

```
i <cr> device information
```

### Example

```
i <cr>
```

### Response

```
?i,pH,2.16 <cr>  
*OK <cr>
```

## Response breakdown

?i,	pH,	2.16
	↑	↑
	Device	Firmware



# Response codes

## Command syntax

**\*OK,1** <cr> enable response **default**

**\*OK,0** <cr> disable response

**\*OK,?** <cr> response on/off?

## Example

## Response

**R** <cr>

**9.560** <cr>  
**\*OK** <cr>

**\*OK,0** <cr>

no response, **\*OK** disabled

**R** <cr>

**9.560** <cr> **\*OK** disabled

**\*OK,?** <cr>

**?\*OK,1** <cr> or **?\*OK,0** <cr>

## Other response codes

**\*ER** unknown command  
**\*OV** over volt ( $VCC \geq 5.5V$ )  
**\*UV** under volt ( $VCC \leq 3.1V$ )  
**\*RS** reset  
**\*RE** boot up complete, ready  
**\*SL** entering sleep mode  
**\*WA** wake up

These response codes  
cannot be disabled

# Reading device status

## Command syntax

Status <cr> voltage at Vcc pin and reason for last restart

### Example

```
Status <cr>
```

### Response

```
?Status,P,5.038 <cr>  
*OK <cr>
```

## Response breakdown

?Status,	P,	5.038
	↑	↑
	Reason for restart	Voltage at Vcc

### Restart codes

P	powered off
S	software reset
B	brown out
W	watchdog
U	unknown

# Sleep mode/low power

## Command syntax

Send any character or command to awaken device.

Sleep <cr> enter sleep mode/low power

## Example

## Response

Sleep <cr>

\*OK <cr>

\*SL <cr>

Any command

\*WA <cr> wakes up device

5V

STANDBY

16 mA


SLEEP

1.16 mA


3.3V

13.9 mA


0.995 mA



Standby  
16 mA



Sleep <cr>



Sleep  
1.16 mA

# Factory reset

## Command syntax

Clears calibration  
LED on  
"\*OK" enabled

Factory <cr> enable factory reset



## Example

## Response


Factory <cr>

\*OK <cr>

Factory <cr>



(reboot)



\*OK <cr>

\*RS <cr>

\*RE <cr>


Baud rate will not change

# Calibration theory


The accuracy of your readings is directly related to the quality of your calibration.  
(Calibration is not difficult, and a little bit of care goes a long way).

## Single, Two point, or Three point calibration accuracy

Single point calibration




Two point calibration



Three point calibration



## Confirm the pH probe is working correctly



A new Atlas Scientific pH probe, still in its soaker bottle will read a pH of **~3.4-3.8**

If your pH probe gives a reading of **zero**, **seven** or **14** continuously and that reading cannot be changed no matter what solution the probe is in, your probe cannot be calibrated and may be damaged.

Contact Atlas Scientific customer support for assistance.


# Calibration order

If this is your first time calibrating the EZO Complete-pH, we recommend following this calibration order.



## Calibration solutions


The Atlas Scientific EZO Complete-pH can work with any brand or value of calibration solution. **We recommend using calibration solutions that have simple values.**



While you can use calibration solutions that have complex values, we recommend avoiding unnecessary complexity. **Unusually specific calibration values should be treated with suspicion.**

# Best practices for calibration

Always watch the readings throughout the calibration process.  
Issue calibration commands once the readings have stabilized.



**⚠ Never do a blind calibration! ⚠**

Issuing a calibration command before the readings stabilize will result in drifting readings.



# Best practices for calibration

**Avoid extended stabilization time.**



Letting the probes pre-calibration readings stabilize over an extended period will cause your calibrated readings to take a long time to stabilize.




**Avoid frequent recalibrations.**

if it ain't broke, don't fix it.

pH probes lose accuracy slowly. Frequent recalibrations to insure high accuracy will often have the opposite effect. It is far more likely that you will miscalibrate the probe rather than improve its accuracy.





# 1. Mid point calibration





✓ Mid point calibrated

# 2. Low point calibration



✓ Low point calibrated

# 3. High point calibration



✓ High point calibrated

## Optional steps:


Confirm your calibration accuracy using the slope command.  
Recalibrate a single point if required.

# Understanding pH slope

The slope function is a powerful tool used to verify calibration and determine the overall health of a pH probe. By evaluating the slope of a pH probe's response curve, you can determine how well a pH probe was calibrated or when that probe is reaching end of life.

**Slope and calibration are directly related. The slope is updated when a calibration command is given. The slope does not update automatically.**

Generally speaking, all pH probes behave the same way. This means a probe's response to calibration can be compared to a simulated pH probe that is mathematically perfect in all ways.



The slope is broken into three sections; acid, base, and neutral. Each section is evaluated separately.

- Acid (pH 1–6.9)
- Base (pH 7.1–14)
- Neutral (pH 7)


An uncalibrated pH probe will have a mathematically perfect slope. Because no pH probe is mathematically perfect, the slope can be used to determine if the pH was calibrated.

**Uncalibrated slope:** 100, 100, 0 (acid, base, neutral)  
% % mV

The first two numbers are percentages, and the third is millivolts. The slope shows that the probe's response to acid and base is 100% correct, and it detects 0 mV in a pH 7. Because such perfection does not exist in the real world, we know this probe was not calibrated.

# Understanding pH slope

pH 7 is the absence of pH; it is not an acid or a base. Therefore it should always be your first calibration point. It is equivalent to the tare function on a scale because it establishes the probe's zero point.



**cal,mid,7**  
**\*OK**

After pH 7 calibration, use the slope command to see how the probe performed during calibration.

**The slope after pH 7 calibration:** 100, 100, -1.2

Here we see the probe reads -1.2mV in pH 7. The closer this number is to 0, the better. A new pH probe should give a millivolt offset no greater than -5mV to 5mV. Over time this number's distance to 0 may increase; the larger the number, the lower the accuracy. A reading >10mV will result in noticeable performance issues.

**It is important to remember that a high number is not definitive evidence that the probe is inaccurate or malfunctioning. It is very common to see a high number if the calibration solution was contaminated and not actually its stated value.**

# Understanding pH slope

The next two calibration points (*pH 4* and *pH 10*) report their slope in percentage. A new pH probe should have a slope of  $>95\%$ .


**The slope after pH 4 calibration:** 98.2, 100, -1.2

**The slope after pH 10 calibration:** 98.2, 97.8, -1.2

## Tips:

Throughout this explanation, we have looked at the slope after each calibration event. This is unnecessary; in reality, it is best to fully calibrate the probe and look at the slope once calibration has been completed.

To gain a deeper understanding of how slope affects the stability and accuracy of a pH probe, intentionally miscalibrate the probe and see how it affects the slope.



# Datasheet change log

## **Datasheet V 1.1**

Revised probe artwork.

## **Datasheet V 1.0**

Revised entire document.

# Firmware updates

V1.5 – Baud rate change (Nov 6, 2014)

- Change default baud rate to 9600

# Warranty

Atlas Scientific™ Warranties the EZO Complete device to be free of defects during the debugging phase of device implementation or 30 days after receiving the EZO Complete device (*whichever comes first*).

## The debugging phase

As defined by Atlas Scientific™, the debugging phase is when the EZO Complete device is connected to a computer to evaluate its output and/or is being integrated into custom software.

**The following activities will void the EZO Complete device warranty:**

- **Soldering any part of the EZO™ class device.**
- **Removing any potting compound.**
- **Embedding the EZO Complete device into a custom machine.**

## Reasoning behind this warranty

**Atlas Scientific™ does not sell consumer electronics.** Once the device has been embedded into a custom-made machine, Atlas Scientific™ cannot possibly warranty the EZO Complete device against the thousands of possible variables that may cause the device to malfunction.

## Please keep this in mind:

- 1. All Atlas Scientific™ devices have been designed to be embedded into a custom-made machine by you, the embedded systems engineer.**
- 2. All Atlas Scientific™ devices have been designed to run indefinitely without failure in the field.**

Atlas Scientific™ is simply stating that once the device is being used in your machine or application, Atlas Scientific™ can no longer take responsibility for the device's continued operation. Doing so would be equivalent to Atlas Scientific™ taking responsibility for the correct operation of your entire machine.